

with promising results in terms of infection eradication and antibiotic resistance selection control.<sup>10</sup> In this scenario, the novel commercial MG + parC (beta) assay, demonstrating to be rapid, simple and accurate, could complement this resistance-guided therapy approach with the additional detection of fluoroquinolone resistance-associated mutations in *M. genitalium*, optimizing and refining antimicrobial stewardship. Nevertheless, despite the test detects known fluoroquinolone resistance-associated SNPs, further studies are required to fully establish the contribution of other newly parC mutations in fluoroquinolone resistance.<sup>6,7</sup> Moreover, more investigations on new antibiotics and novel combinations with existing treatments are imperative to fight against *M. genitalium* that may soon become the next sexually transmitted superbug.

## Funding

SpeeDx Pty Ltd. supplied all the reagents for molecular testing of *M. genitalium* and fluoroquinolone resistance.

## Conflicts of interest

J.S.P and J.E. employed by the "Vall d'Hebron University Hospital", have received remuneration for contract work from SpeeDx Pty Ltd. M.E. employed by the "Corporació Sanitària Parc Taulí", has received remuneration for contract work from SpeeDx Pty Ltd. M.F.H. and M.E. have participated in symposiums organized by SpeeDx Pty Ltd.

## References

- Taylor-Robinson D, Jensen JS. *Mycoplasma genitalium*: from Chrysalis to multi-colored butterfly. *Clin Microbiol Rev*. 2011;24:498–514.
- Jensen JS, Cusini M, Gomberg M, Moi H. 2016 European guideline on *Mycoplasma genitalium* infections. *J Eur Acad Dermatol Venereol*. 2016;30:1686–93.
- Lau A, Bradshaw CS, Lewis D, Fairley CK, Chen MY, Kong FY, et al. The efficacy of azithromycin for the treatment of genital *Mycoplasma genitalium*: a systematic review and meta-analysis. *Clin Infect Dis*. 2015;61:1389–99.

## Application of a new method for the detection of catheter colonization and catheter-related bacteraemia in newborns

### Aplicación de un nuevo método para la detección de colonización y bacteriemia relacionada con catéter en neonatos

Healthcare associated infections (HAIs) in neonatal units are increasing due to advances in invasive therapeutic and diagnostic procedures and increased survival of preterm babies being bloodstream infections (BSIs) one of the most common type of HAIs.<sup>1,2</sup> The drawing of paired blood samples from the catheter and a peripheral vein is considered the best approach to reduce suspicion of colonization of the device, contamination of the sample or true bloodstream infection caused by common commensal microorganisms.<sup>3</sup>

Current guidelines for the diagnosis of catheter-related bloodstream infections (CRBSI) recommend culture of catheters only when clinical suspicion are present to avoid unnecessary antimicrobial treatment.<sup>4</sup> Catheter colonization can be diagnosed either by rolling the catheter tip on the agar surface (Maki's roll-plate technique, MT) or by sonicating the distal 4–5 cm of the catheter.<sup>5</sup> However, MT is not fully reliable in the case of silicone neonatal peripherally inserted central catheter (SN-PICC) tips based on the study of Martín-Rabandan et al. due to the low capacity of this method for microbial detection inside catheters. Culture results may improve after catheter fragmentation.<sup>6</sup> In fact, the new protocol of 'Microbiological diagnosis of infections associated with

- Barberá MJ, Fernández-Huerta M, Jensen JS, Caballero E, Andreu A. *Mycoplasma genitalium* macrolide and fluoroquinolone resistance: prevalence and risk factors among a 2013–2014 cohort of patients in Barcelona, Spain. *Sex Transm Dis*. 2017;44:457–62.
- Piñeiro L, Idigoras P, de la Caba I, López-Olaizola M, Cilla G. Guided antibiotic therapy for *Mycoplasma genitalium* infections: analysis of mutations associated with resistance to macrolides and fluoroquinolones. *Enferm Infect Microbiol Clin*. 2019;37:394–7.
- Murray GL, Bradshaw CS, Bissessor M, Danielewski J, Twin J, Jensen JS, et al. Increasing macrolide and fluoroquinolone resistance in *Mycoplasma genitalium*. *Emerg Infect Dis*. 2017;23:809–12.
- Couldwell DL, Tagg KA, Jeoffreys NJ, Gilbert GL. Failure of moxifloxacin treatment in *Mycoplasma genitalium* infections due to macrolide and fluoroquinolone resistance. *Int J STD AIDS*. 2013;24:822–8.
- Fernández-Huerta M, Bodiyabadu K, Esperalba J, Bradshaw CS, Serra-Pladevall J, Garland SM, et al. Multicenter clinical evaluation of a novel multiplex real-time (qPCR) assay for detection of fluoroquinolone resistance in *Mycoplasma genitalium*. *J Clin Microb*. 2019 [Epub ahead of print].
- Fernández-Huerta M, Vall M, Fernández-Naval C, Barberá MJ, Arando M, López L, et al. *Mycoplasma genitalium* macrolide resistance update: rate among a 2016–2017 cohort of patients in Barcelona, Spain. *Enferm Infect Microbiol Clin*. 2019 [Epub ahead of print].
- Read TRH, Fairley CK, Murray GL, Jensen JS, Danielewski J, Worthington K, et al. Outcomes of resistance-guided sequential treatment of *Mycoplasma genitalium* infections: a prospective evaluation. *Clin Infect Dis*. 2018;68:554–60.

Miguel Fernández-Huerta <sup>a</sup>, Juliana Esperalba <sup>a</sup>,  
Judit Serra-Pladevall <sup>a</sup>, Mateu Espasa <sup>b,\*</sup>

<sup>a</sup> Microbiology Department, Vall d'Hebron University Hospital, Universitat Autònoma de Barcelona, Barcelona, Spain

<sup>b</sup> Microbiology Department, Corporació Sanitària Parc Taulí, Universitat Autònoma de Barcelona, Barcelona, Spain

\* Corresponding author.

E-mail address: mespasa@tauli.cat (M. Espasa).

<https://doi.org/10.1016/j.eimc.2019.10.003>

0213-005X/ © 2019 Elsevier España, S.L.U. and Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica. All rights reserved.

intravascular catheters' elaborated by the Spanish Society of Infectious Diseases and Clinical Microbiology (SEIMC) includes roll-plate after slicing (RPS) as a good option in CRBSI diagnosis.<sup>7</sup>

Our objective was to compare the yield of MT and RPS for the detection of colonization and CRBSI. For this purpose, tips of neonatal catheters, all of them SN-PICCs, were prospectively cultured in the Hospital Miguel Servet of Zaragoza for 2 years (January 1, 2017 to December 31, 2018). Once the MT was performed, the catheter was opened longitudinally and the contents were extended with swab, leaving the catheter in the culture plate. They were incubated 4 days, issuing a preliminary report at 48 h. In all cases, culture was carried out for *Malassezia furfur*, adding a drop of olive oil in the initial area of the seed preventing it from spreading.

It was considered that the catheter was colonized when either of the two methods found ≥15 colony forming units (CFU). CRBSI was defined when the same microorganism was isolated in the catheter culture and peripheral blood cultures, without another source of infection.

During the study period, 79 sPICCS were processed. The number of colonized catheters was 24 (30.4%), with TM detected 15/24 and with RPS 21/24 ( $p = 0.146$ ). The microorganisms isolated in colonized catheters were 18 Gram positive, 4 Gram-negative and 7 yeast, 5 of them *M. furfur*. BRC was diagnosed in 11 patients (13.9%), in 6 of them the diagnosis was made only with RPS ( $p = 0.031$ ). Isolated microorganisms are shown in Table 1, in which the CFUs are expressed as: negative, >15 when the CFUs are between 15–100 and >100.

**Table 1**

Isolated microorganisms using MT and RPS.

	Maki's roll-plate technique	CFU	Roll-plate after slicing	CFU	CRBSI
1	<i>S. epidermidis</i>	>100	<i>S. epidermidis</i>	>100	Yes
2	<i>S. epidermidis</i> + <i>S. haemolyticus</i>	>15/>100	<i>S. epidermidis</i> + <i>S. haemolyticus</i>	>15/>100	Yes
3	Negative		<i>S. epidermidis</i>	>15	Yes
4	Negative		<i>S. epidermidis</i>	>100	Yes
5	<i>S. epidermidis</i>	<15	<i>S. epidermidis</i>	>15	Yes
6	<i>S. haemolyticus</i>	15	<i>S. haemolyticus</i>	>15	Yes
7	Negative		<i>S. aureus</i>	15	Yes
8	<i>E. faecalis</i>	<15	<i>E. faecalis</i>	>100	Yes
9	<i>E. cloacae</i> complex	>15	<i>E. cloacae</i> complex	>100	Yes
10	<i>E. cloacae</i> complex + <i>K. oxytoca</i>	<15/<15	<i>E. cloacae</i> complex + <i>K. oxytoca</i>	>100/>15	Yes
11	<i>C. lusitaniae</i> + <i>C. lypolitica</i>	>100/>15	<i>C. lusitaniae</i> + <i>C. lypolitica</i>	>100/>15	Yes
12	<i>S. epidermidis</i>	>15	Negative		No
13	<i>S. haemolyticus</i>	>15	Negative		No
14	<i>S. epidermidis</i>	>15	<i>S. epidermidis</i>	<15	No
15	Negative		<i>S. hominis</i>	>15	No
16	<i>S. epidermidis</i>	>15	<i>S. epidermidis</i>	>15	No
17	<i>S. epidermidis</i>	>100	<i>S. epidermidis</i>	>15	No
18	<i>S. epidermidis</i>	>15	<i>S. epidermidis</i>	>15	No
19	<i>S. epidermidis</i>	>100	<i>S. epidermidis</i>	>100	No
20	Negative		<i>M. furfur</i> + <i>S. marcescens</i>	>100/<15	No
21	Negative		<i>M. furfur</i>	>100	No
22	<i>M. furfur</i>	>15	<i>M. furfur</i>	>100	No
23	<i>M. furfur</i>	>100	<i>M. furfur</i>	>100	No
24	<i>S. epidermidis</i>	>15	<i>S. epidermidis</i> + <i>M. furfur</i>	>100/>15	No



**Fig. 1.** On the left, a catheter seeded by MT with <15 CFU of *E. cloacae* complex and *K. oxytoca*. On the right, RPS with significant count of *E. cloacae* complex (>15 CFU) and *K. oxytoca* (>100 UFC).

In our hospital, CRBSI is an important cause of nosocomial infection. For the diagnosis of colonization and CRBSI, RPS presented significantly higher returns than MT, so the data suggests that the diagnosis of CRBSI would be improved with this new technique being interesting to use RPS routinely. Currently, in our laboratory the two techniques are being carried out in parallel to increase the sample size and obtain stronger conclusions (Fig. 1).

### Conflict of interest

The authors have no conflict of interest or source of funding.

### References

- Kumar S, Shankar B, Arya S, Deb M, Chellani H. Healthcare associated infections in neonatal intensive care unit and its correlation with environmental surveillance. *J Infect Public Health*. 2018;11:275–9.
- Cho HJ, Cho HK. Central line-associated bloodstream infections in neonates. *Korean J Pediatr*. 2019;62:79–84.
- Ferreira J, Camargos PAM, Clemente WT, Romanelli RMC. Clinical usefulness of catheter-drawn blood samples and catheter tip cultures for the diagnosis of catheter-related bloodstream infections in neonatology: a systematic review. *Am J Infect Control*. 2018;46:81–7.
- Kitano T, Takagi K, Arai I, Yasuhara H, Ebisu R, Ohgitani A, et al. Efficacy of routine catheter tip culture in a neonatal intensive care unit. *Pediatr Int*. 2018;60:423–7.
- Guembe M, Martín-Rabadán P, Cruces R, Pérez Granda MJ, Bouza E. Slicing silicone neonatal vascular catheter tips improves colonization detection by the roll-plate technique. *Clin Microbiol Infect*. 2017;23:410, e1–410.e3.
- Martín-Rabadán P, Pérez García F, Zamora Flores E, Nisa ES, Guembe M, Bouza E. Improved method for the detection of catheter colonization and catheter-related bacteremia in newborns. *Diagn Microbiol Infect Dis*. 2017;87:311–4.
- Aldea C, Gracia I, Guembe M, Martínez J. Diagnóstico microbiológico de las infecciones asociadas a catéteres intravasculares. 15a. Procedimientos en Microbiología Clínica. Editores: Cercenado E, Cantón R. SEIMC. 2018.

Saray Mormeneo Bayo\*, María Pilar Palacián Ruíz,  
Diego Ortega Larrea, María Cruz Villuendas Usón

Servicio de Microbiología, Hospital Universitario Miguel-Servet,  
Zaragoza, Instituto de Investigación Sanitaria Aragón, Spain

\* Corresponding author.

E-mail address: samorbayo@hotmail.com (S. Mormeneo Bayo).

<https://doi.org/10.1016/j.eimc.2019.10.004>

0213-005X/ © 2019 Elsevier España, S.L.U. and Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica. All rights reserved.